

## Product datasheet for **MP204358**

### Etv1 Mouse qPCR Primer Pair (NM\_007960)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14009
Forward Sequence:	TCCTGGCTCATCCAAGCAGAAC
Reverse Sequence:	CGGTACATTCCAGGCTCTTGCT
ACCN:	<u><a href="#">BC005645</a></u> , <u><a href="#">NM_007960</a></u> , <u><a href="#">NM_007960.1</a></u> , <u><a href="#">NM_007960.2</a></u> , <u><a href="#">NM_007960.3</a></u> , <u><a href="#">NM_007960.4</a></u> , <u><a href="#">NM_007960.5</a></u> , <u><a href="#">BB751428</a></u> , <u><a href="#">BE859344</a></u> , <u><a href="#">BY062397</a></u> , <u><a href="#">BY122680</a></u> , <u><a href="#">BY330680</a></u>
UniProt ID:	<u><a href="#">P41164</a></u>
Synonyms:	ER81; Etsrp81
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)