

Product datasheet for **MP204390**

Eif2b1 Mouse qPCR Primer Pair (NM_145371)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	209354
Forward Sequence:	GAGTGGTAGAGAACGGAGGCAT
Reverse Sequence:	CTCCTGGTTGAGTGGAAGAGG
ACCN:	<u>BC003426, NM_145371, NM_145371.1, NM_145371.2, NM_145371.3, NM_145371.4</u>
UniProt ID:	<u>Q99LC8</u>
Synonyms:	26kD; D5Ertd406; D5Ertd406e; EIF2; EIF2BA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)