

Product datasheet for **MP204429**

Eif5a Mouse qPCR Primer Pair (NM_181582)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	276770
Forward Sequence:	AGAGGACCTTCGTCTGCCTGAA
Reverse Sequence:	TGTCATGGCAGACAGCACTGTG
ACCN:	<u>BC003889</u> , <u>BC008093</u> , <u>BC024899</u> , <u>NM_181582</u> , <u>NM_181582.1</u> , <u>NM_181582.2</u> , <u>NM_181582.3</u> , <u>NM_181582.4</u>
UniProt ID:	<u>P63242</u>
Synonyms:	AA410058; D19Wsu54; D19Wsu54e; eIF-4D; eIF-5A; eIF-5A-1; eIF-5A1; Eif4d; Eif5a1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)