

## Product datasheet for **MP204460**

### Exosc1 Mouse qPCR Primer Pair (NM\_025644)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66583
Forward Sequence:	GCCTGATGAAGACCAGCGAGAA
Reverse Sequence:	GGCAAACCGTGAGTTGATGCTAG
ACCN:	<u><a href="#">BC024423</a>, <a href="#">NM_025644</a>, <a href="#">NM_025644.1</a>, <a href="#">NM_025644.2</a>, <a href="#">NM_025644.3</a>, <a href="#">NM_025644.4</a>, <a href="#">BY361617</a></u>
UniProt ID:	<u><a href="#">Q9DAA6</a></u>
Synonyms:	2610035C18Rik; 2610104C07Rik; 2610312F07Rik; AI447561
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)