

Product datasheet for **MP204461**

Exosc10 Mouse qPCR Primer Pair (NM_016699)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	50912
Forward Sequence:	CTGTCCAGTTGCAGGTGGTTTG
Reverse Sequence:	TGCTGCGAGTTCAGGTGCTTCT
ACCN:	<u>BC027326</u> , <u>BC027840</u> , <u>NM_016699</u> , <u>NM_016699.1</u> , <u>NM_016699.2</u> , <u>BC003952</u> , <u>BY277385</u> , <u>NM_016699.3</u>
UniProt ID:	<u>P56960</u>
Synonyms:	p2; p3; p4; PM-Scl; PM/Scl-100; Pmscl2; RRP6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)