

Product datasheet for **MP204561**

Tafa5 Mouse qPCR Primer Pair (NM_134096)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	106014
Forward Sequence:	GCTCAAGTCTGTCCAGCTGAAC
Reverse Sequence:	AAGCGGCTAGTGAGGATGGATG
ACCN:	BC015306 , NM_134096 , NM_134096.1 , NM_134096.2 , BB873716 , BE945208 , BI986174 , BU057253 , BY733542
UniProt ID:	Q91WE9
Synonyms:	Fam19; Fam19a5; Tafa-5; Tar; Tara-5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)