

Product datasheet for **MP204614**

Fancd2 Mouse qPCR Primer Pair (NM_001033244)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 211651 |
| Forward Sequence: | TCAACCTCCTGCCGCTGTTCTA |
| Reverse Sequence: | TTGCCACGCAAAGTCTCAGCAG |
| ACCN: | <u>NM_001033244</u> , <u>NM_001033244.1</u> , <u>NM_001033244.2</u> , <u>NM_001033244.3</u> , <u>BC139004</u> , <u>BC040776</u> , <u>BC042619</u> , <u>BC042870</u> |
| UniProt ID: | <u>Q80V62</u> |
| Synonyms: | 2410150O07Rik; AU015151; BB137857; FA-D2; FA4; FACD; FAD; FANCD |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)