

Product datasheet for **MP204632**

FasL Mouse qPCR Primer Pair (NM_010177)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14103
Forward Sequence:	GAAGGAACTGGCAGAACTCCGT
Reverse Sequence:	GCCACACTCCTCGGCTCTTTT
ACCN:	<u>BC052866</u> , <u>NM_010177</u> , <u>NM_010177.1</u> , <u>NM_010177.2</u> , <u>NM_010177.3</u> , <u>NM_010177.4</u>
UniProt ID:	<u>P41047</u>
Synonyms:	APT1LG1; CD95-L; CD95L; CD178; Fas-L; Faslg; gld; Tnfsf6; Tnlg1a
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)