

Product datasheet for **MP204653**

Fabp3 Mouse qPCR Primer Pair (NM_010174)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14077
Forward Sequence:	AGAGTTCGACGAGGTGACAGCA
Reverse Sequence:	TTGTCTCCTGCCCGTTCCACTT
ACCN:	<u>BC002082</u> , <u>BC089542</u> , <u>NM_010174</u> , <u>NM_010174.1</u> , <u>NM_010174.2</u>
UniProt ID:	<u>P11404</u>
Synonyms:	Fabph-1; Fabph-4; Fabph1; Fabph4; H-FABP; Mdgi
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)