

Product datasheet for **MP204678**

Fbxl3 Mouse qPCR Primer Pair (NM_015822)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	50789
Forward Sequence:	GAACTGAATCAGCCAGCGACATC
Reverse Sequence:	GCGACAGTATGTACAAGCAGC
ACCN:	BC067203 , NM_015822 , NM_015822.1 , NM_015822.2 , BC026835 , BC037665 , BC046330 , BY146960 , BY276318 , NM_015822.3
UniProt ID:	Q8C4V4
Synonyms:	AU041772; AW212966; FBK; Fbl3a; Fbxl3a; Ovtm; Play68
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)