

## Product datasheet for **MP204684**

### Fbxw10 Mouse qPCR Primer Pair (NM\_001033669)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	213980
Forward Sequence:	GGGACAATCACTTGCCTGGATG
Reverse Sequence:	AGGTCTCACTGATCTTGGCAGC
ACCN:	<u><a href="#">NM_001033669</a></u> , <u><a href="#">NM_001033669.1</a></u> , <u><a href="#">NM_001033669.2</a></u> , <u><a href="#">BC085108</a></u> , <u><a href="#">BC141212</a></u> , <u><a href="#">BC172760</a></u> , <u><a href="#">BX522544</a></u> , <u><a href="#">BY360237</a></u>
Synonyms:	Fbw10; HREP; SM2SH2; SM25H2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)