

Product datasheet for **MP204699**

Fbxo11 Mouse qPCR Primer Pair (NM_001081034)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	225055
Forward Sequence:	TCAAGCAGGTGTCCTCATCAGC
Reverse Sequence:	CTGATTGCCTTCGAGTGTTCAG
ACCN:	<u>BC117884</u> , <u>BC117885</u> , <u>BC128479</u> , <u>NM_001081034</u> , <u>NM_001081034.1</u> , <u>NM_001081034.2</u> , <u>BC027802</u> , <u>BC049946</u> , <u>BC055343</u> , <u>BC110686</u> , <u>NM_001081034.3</u>
UniProt ID:	<u>Q7TPD1</u>
Synonyms:	C80048; Jf
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)