

Product datasheet for **MP204713**

Fbxo32 Mouse qPCR Primer Pair (NM_026346)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67731
Forward Sequence:	CTTCTCGACTGCCATCCTGGAT
Reverse Sequence:	TCTTTTGGGCGATGCCACTCAG
ACCN:	<u>BC027211</u> , <u>NM_026346</u> , <u>NM_026346.1</u> , <u>NM_026346.2</u> , <u>BB838158</u> , <u>BU921446</u>
UniProt ID:	<u>Q9CPU7</u>
Synonyms:	4833442G10Rik; AI430017; ATROGIN1; Gm20361; MAFbx
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)