

Product datasheet for **MP204720**

Fbxo41 Mouse qPCR Primer Pair (NM_001001160)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	330369
Forward Sequence:	TGGCTTCACTCAAGCAGGACCT
Reverse Sequence:	CTCCTCGATGAAGTGTGCAGC
ACCN:	<u>BC070445</u> , <u>NM_001001160</u> , <u>NM_001001160.1</u> , <u>NM_001001160.2</u> , <u>NM_001001160.3</u> , <u>BC048554</u> , <u>BC061475</u>
UniProt ID:	<u>Q6NS60</u>
Synonyms:	9630017H13; D6Ertd538e
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)