

Product datasheet for **MP204776**

Ftl1 Mouse qPCR Primer Pair (NM_010240)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14325
Forward Sequence:	CCTCGAGTTTCAGAACGATCGC
Reverse Sequence:	CCTGATTCAGGTTCTTCCATG
ACCN:	<u>BC019840</u> , <u>BC081462</u> , <u>BC083350</u> , <u>BC085309</u> , <u>BC092259</u> , <u>BC106145</u> , <u>BC106146</u> , <u>NM_010240</u> , <u>NM_010240.1</u> , <u>NM_010240.2</u> , <u>BC020403</u> , <u>BY134899</u>
Synonyms:	Ftl; Ftl-1; L-ferritin
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)