

Product datasheet for **MP204804**

Fgf12 Mouse qPCR Primer Pair (NM_183064)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14167
Forward Sequence:	GGAGAAGGCTATCTCTACAGCTC
Reverse Sequence:	CTGATTCTGCTGGCGATACAG
ACCN:	<u>BC030485</u> , <u>NM_183064</u> , <u>NM_183064.1</u> , <u>NM_183064.2</u> , <u>NM_183064.3</u> , <u>NM_183064.5</u> , <u>BB313752</u> , <u>BB611389</u> , <u>BM941461</u> , <u>BP768136</u> , <u>BY477982</u>
UniProt ID:	<u>P61329</u>
Synonyms:	AV114868; B230343J05Rik; FGF-12; FHF-1; Fhf1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)