

Product datasheet for **MP204865**

Flna Mouse qPCR Primer Pair (NM_010227)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	192176
Forward Sequence:	CAGCAAGCTACAGGTGGAACCT
Reverse Sequence:	TCAGTGGTTCCTCTCGGAAGA
ACCN:	<u>NM_010227</u> , <u>NM_010227.1</u> , <u>NM_010227.2</u> , <u>NM_010227.3</u> , <u>BC151024</u> , <u>BC004061</u> , <u>BC038478</u> , <u>BC039208</u> , <u>BC054432</u> , <u>BC099506</u> , <u>BQ950333</u>
Synonyms:	ABP-280; Dilp2; F730004A14Rik; filamin-1; Fln1; GENA 379
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)