

Product datasheet for **MP204932**

Foxg1 Mouse qPCR Primer Pair (NM_008241)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	15228
Forward Sequence:	TACTACCGCGAGAACAAGCAGG
Reverse Sequence:	GAGCATCCAGTAGTTGCCCTTG
ACCN:	<u>BC046958</u> , <u>NM_008241</u> , <u>NM_008241.1</u> , <u>NM_008241.2</u> , <u>BC055774</u> , <u>BC064449</u> , <u>BC079597</u>
UniProt ID:	<u>Q60987</u>
Synonyms:	2900064B05Rik; BF-1; Bf1; Hfh9; Hfhbf1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)