

## Product datasheet for **MP204952**

### Foxp1 Mouse qPCR Primer Pair (NM\_053202)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	108655
Forward Sequence:	CATGCCTCTACCAATGGACAGC
Reverse Sequence:	GAAGTCGTCACAAACCGCCTCA
ACCN:	<a href="#">BC064764</a> , <a href="#">NM_053202</a> , <a href="#">NM_053202.1</a> , <a href="#">NM_053202.2</a> , <a href="#">BC065386</a>
UniProt ID:	<a href="#">P58462</a>
Synonyms:	3110052D19Rik; 4932443N09Rik; AI461938; AW494214
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)