

## Product datasheet for **MP205021**

### Fuca1 Mouse qPCR Primer Pair (NM\_024243)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	71665
Forward Sequence:	AAGCCCTGTGTCTTGGAACTGG
Reverse Sequence:	TTCCAAGAGCGAGTGGTAGAGG
ACCN:	<u><a href="#">BC003235</a>, <a href="#">NM_024243</a>, <a href="#">NM_024243.1</a>, <a href="#">NM_024243.2</a>, <a href="#">NM_024243.3</a>, <a href="#">NM_024243.4</a>, <a href="#">BY158254</a></u>
UniProt ID:	<u><a href="#">Q99LJ1</a></u>
Synonyms:	0610006A03Rik; 9530055J05Rik; Afuc; Fuca
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)