

Product datasheet for **MP205079**

Hcar2 Mouse qPCR Primer Pair (NM_030701)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80885
Forward Sequence:	CTGTTTCCACCTCAAGTCCTGG
Reverse Sequence:	CATAGTTGTCCGTCAGGAACGG
ACCN:	<u>BC120715</u> , <u>BC125296</u> , <u>NM_030701</u> , <u>NM_030701.1</u> , <u>NM_030701.2</u> , <u>NM_030701.3</u> , <u>BY145731</u> , <u>BY516126</u>
UniProt ID:	<u>Q9EP66</u>
Synonyms:	Gpr109a; Gpr109b; HM74; mHM74b; Niacr1; PUMA-G; Pumag
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)