

Product datasheet for **MP205136**

Gpr37l1 Mouse qPCR Primer Pair (NM_134438)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	171469
Forward Sequence:	GATTACCAAGCAAAGGCTACTCG
Reverse Sequence:	GCCACATGGAATCGGTCTATGC
ACCN:	BC098482 , BC127619 , BC127620 , NM_134438 , NM_134438.1 , NM_134438.2 , NM_134438.3 , BY727284
UniProt ID:	Q99JG2
Synonyms:	AW047233; CAG-18; D0Kist8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)