

Product datasheet for **MP205179**

Gabpb2 Mouse qPCR Primer Pair (NM_029885)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	213054
Forward Sequence:	GACCAAGGATGGCAGAGATGAC
Reverse Sequence:	CTTGAGGAGCTGGTGTCTGTAC
ACCN:	<u>NM_029885</u> , <u>NM_029885.1</u> , <u>BC150160</u> , <u>BC150482</u> , <u>BP768545</u>
UniProt ID:	<u>P81069</u>
Synonyms:	1810015F01Rik; 5830427M07; 9430006E19Rik; A430024B14Rik; AV050852; Gabpb2-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)