

Product datasheet for **MP205184**

Galt Mouse qPCR Primer Pair (NM_016658)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14430
Forward Sequence:	GCAGACCTATCACAGCCAGCAT
Reverse Sequence:	TCCAGTGCTCACTGGTTAGGAC
ACCN:	BC010985 , NM_016658 , NM_016658.1 , NM_016658.2 , NM_016658.3 , BB610075 , BB868849 , BY282576 , BY289084 , BY416479
UniProt ID:	Q03249
Synonyms:	AW553376
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)