

Product datasheet for **MP205231**

Ggt5 Mouse qPCR Primer Pair (NM_011820)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23887
Forward Sequence:	AACACCCTGGAGACTGTAGCCA
Reverse Sequence:	AAGGCAGCCAAGTCCTGAACTG
ACCN:	BC113775 , NM_011820 , NM_011820.1 , NM_011820.2 , NM_011820.3 , NM_011820.4 , BC051971 , BC089585 , NM_011820.5
UniProt ID:	Q9Z2A9
Synonyms:	AI551243; GGL; GGT-REL; Ggtla1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)