

Product datasheet for **MP205232**

Ggt6 Mouse qPCR Primer Pair (NM_027819)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	71522
Forward Sequence:	GCAATGTCGTGGATGCTGGAGT
Reverse Sequence:	GGCAGTTGAGTTTCCTGAAGAGC
ACCN:	<u>BC055837</u> , <u>BC058747</u> , <u>NM_027819</u> , <u>NM_027819.1</u> , <u>NM_027819.2</u>
UniProt ID:	<u>Q6PDE7</u>
Synonyms:	9030405D14Rik; AI427679; AW490214
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)