

Product datasheet for **MP205285**

Gpbp1 Mouse qPCR Primer Pair (NM_028487)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	73274
Forward Sequence:	GAGAACGCAGACAGTTTGAGGC
Reverse Sequence:	GGTTGGGTATGTGTGTCTGG
ACCN:	BC016083 , BC067075 , BC094888 , NM_028487 , NM_028487.1 , NM_028487.2 , NM_028487.3 , NM_028487.4 , BC067075.1 , BC076635 , BC096442 , NM_028487.5
UniProt ID:	Q6NXH3
Synonyms:	1700034P14Rik; AU019836; D230035M11Rik; Gpbp; mGPBP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)