

Product datasheet for **MP205426**

Gtf3c2 Mouse qPCR Primer Pair (NM_027901)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	71752
Forward Sequence:	GTAGCCTGGAAGTGCCAGTTGT
Reverse Sequence:	CGGTAGATTGTGCCATCCTCAG
ACCN:	<u>BC034369</u> , <u>BC043100</u> , <u>NM_027901</u> , <u>NM_027901.1</u> , <u>NM_027901.2</u> , <u>BC026573</u>
UniProt ID:	<u>Q8BL74</u>
Synonyms:	1300004C11Rik; 2610510G03Rik; AI225816; AU041069; mKIAA0011; TFIIC110
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)