

Product datasheet for **MP205428**

Gtf3c4 Mouse qPCR Primer Pair (NM_172977)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	269252
Forward Sequence:	GCTGGAGCTTATCTGCGATGTG
Reverse Sequence:	GCGAACTTTTCCTTGCACTCAGC
ACCN:	<u>BC095998</u> , <u>BC117919</u> , <u>BC117920</u> , <u>NM_172977</u> , <u>NM_172977.1</u> , <u>NM_172977.2</u> , <u>NM_172977.3</u> , <u>BC061476</u> , <u>BU610856</u>
UniProt ID:	<u>Q8BMQ2</u>
Synonyms:	5330400C03; AI426938; AU014771; AU017413; KAT12; TFIIC90
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)