

Product datasheet for **MP205553**

GIs Mouse qPCR Primer Pair (NM_001081081)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14660
Forward Sequence:	CAGAAGGCACAGACATGGTTGG
Reverse Sequence:	CAAGGTGGCAGCCATCACA
ACCN:	NM_001081081 , NM_001081081.1 , NM_001081081.2 , BC025534 , BC050208 , BC169203 , BE849049
UniProt ID:	D3Z7P3
Synonyms:	6330442B14; AI314027; B230365M23Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)