

Product datasheet for **MP205565**

Eprs Mouse qPCR Primer Pair (NM_029735)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	107508
Forward Sequence:	GAGCGGAGAAACCAGAAGCCAA
Reverse Sequence:	CGATTCTGTCAGCGACTTGAC
ACCN:	<u>NM_029735</u> , <u>NM_029735.1</u> , <u>BC141049</u> , <u>BC003226</u> , <u>BC040802</u> , <u>BC058337</u> , <u>BC094679</u> , <u>BU704580</u> , <u>NM_029735.2</u>
UniProt ID:	<u>Q8CGC7</u>
Synonyms:	2410081F06Rik; 3010002K18Rik; C79379; Qprs
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)