

Product datasheet for **MP205607**

Glyctk Mouse qPCR Primer Pair (NM_174846)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	235582
Forward Sequence:	ACGATCCAGGAGCTGAACACCA
Reverse Sequence:	CTGGCAATCACTTCCAAAGGGTC
ACCN:	<u>BC025834</u> , <u>BC025935</u> , <u>BC030732</u> , <u>BC033063</u> , <u>BC036136</u> , <u>NM_174846</u> , <u>NM_174846.1</u> , <u>NM_174846.2</u> , <u>NM_174846.3</u> , <u>NM_174846.4</u> , <u>BC036136.1</u> , <u>BB747408</u> , <u>BC030736</u> , <u>BY111608</u>
UniProt ID:	<u>Q8QZY2</u>
Synonyms:	6230410P16Rik; C430006L12
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)