

## Product datasheet for **MP205689**

### Gga3 Mouse qPCR Primer Pair (NM\_173048)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	260302
Forward Sequence:	GACAACCTGTCCAGGGTCATCA
Reverse Sequence:	CGCCTTGGTTACCACAGTGACT
ACCN:	<a href="#">NM_173048</a> , <a href="#">NM_173048.1</a> , <a href="#">NM_173048.2</a> , <a href="#">NM_173048.3</a> , <a href="#">BC020366</a> , <a href="#">BC024266</a> , <a href="#">BC025051</a> , <a href="#">BC056424</a> , <a href="#">BC085617</a> , <a href="#">BC110676</a> , <a href="#">BC148430</a> , <a href="#">BY268727</a> , <a href="#">BY296723</a>
UniProt ID:	<a href="#">Q8BMI3</a>
Synonyms:	C230037M19Rik; mKIAA0154
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)