

Product datasheet for **MP205697**

Golim4 Mouse qPCR Primer Pair (NM_175193)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	73124
Forward Sequence:	GAGCATCTGGAGGAAGAACACG
Reverse Sequence:	TGTAGCTGCCTTGGTCGAATGC
ACCN:	<u>BC086642</u> , <u>BC089354</u> , <u>NM_175193</u> , <u>NM_175193.1</u> , <u>NM_175193.2</u> , <u>NM_175193.3</u> , <u>NM_175193.4</u> , <u>NM_175193.5</u> , <u>NM_175193.6</u> , <u>BC023824</u> , <u>BC049244</u>
UniProt ID:	<u>Q8BXA1</u>
Synonyms:	3110027H23Rik; GIMPC; Golph4; GPP130; P138
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

