

Product datasheet for **MP205699**

Golph3 Mouse qPCR Primer Pair (NM_025673)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66629
Forward Sequence:	AGACTCAGCCTCCAGAGACAGT
Reverse Sequence:	GTTTCTCTGTCGTCAGTACACCC
ACCN:	<u>BC031445</u> , <u>BC051035</u> , <u>NM_025673</u> , <u>NM_025673.1</u> , <u>NM_025673.2</u> , <u>BC031445.1</u> , <u>BC003961</u> , <u>BC027333</u> , <u>BC039909</u> , <u>BC058233</u>
UniProt ID:	<u>Q9CRA5</u>
Synonyms:	4733401N08Rik; 5730410D03Rik; AW413496
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

