

Product datasheet for **MP205702**

Gorasp2 Mouse qPCR Primer Pair (NM_027352)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	70231
Forward Sequence:	TCGGTCACACCGAGCAACTTGT
Reverse Sequence:	CACTGTGAGGTCTCAATCCTGC
ACCN:	<u>BC005600</u> , <u>BC016455</u> , <u>NM_027352</u> , <u>NM_027352.1</u> , <u>NM_027352.2</u> , <u>NM_027352.3</u> , <u>NM_027352.4</u> , <u>BC007468</u> , <u>BC024982</u> , <u>BU707592</u>
UniProt ID:	<u>Q99JX3</u>
Synonyms:	0610011A07Rik; 2410043M02Rik; 5730520M13Rik; 9430094F20Rik; AW552058
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

