

Product datasheet for **MP205725**

Grhl2 Mouse qPCR Primer Pair (NM_026496)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	252973
Forward Sequence:	GGATGTGAACGAGGAGGCAAAG
Reverse Sequence:	CTGGCAGTATGCTCTGTGGATG
ACCN:	BC055035 , NM_026496 , NM_026496.2 , NM_026496.3 , NM_026496.4 , BC004783 , BC020320 , BC041795 , BC042575 , BY018142 , BY703140
UniProt ID:	Q8K5C0
Synonyms:	0610015A08Rik; BOM; clft3; Tcfcp2l3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)