

Product datasheet for **MP205734**

Ccser2 Mouse qPCR Primer Pair (NM_027045)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	72972
Forward Sequence:	GAGTCGTGTGACCTTGAGGATG
Reverse Sequence:	AATCCGTCCTGCGACTGTCGTA
ACCN:	<u>NM_027045</u> , <u>NM_027045.1</u> , <u>BC139090</u> , <u>BB871809</u> , <u>BC020073</u> , <u>BC048641</u> , <u>BC060947</u> , <u>BC111465</u> , <u>BM214039</u>
UniProt ID:	<u>Q3UHI0</u>
Synonyms:	I700012P13Rik; 2900054P12Rik; AI035535; AW108503; Fam190b; Gcap14
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

