

Product datasheet for **MP205747**

Gigyf2 Mouse qPCR Primer Pair (NM_146112)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	227331
Forward Sequence:	AGATGCAGCACCTCCAGTACCA
Reverse Sequence:	TCTTCAGAGCCTGGAAGTCTG
ACCN:	<u>BC089036</u> , <u>NM_146112</u> , <u>NM_146112.1</u> , <u>NM_146112.2</u> , <u>NM_146112.3</u> , <u>NM_146112.4</u> , <u>BC006835</u> , <u>BC021370</u> , <u>BC025466</u> , <u>BC027137</u> , <u>BC030845</u> , <u>BC076610</u> , <u>BC082811</u> , <u>BY117987</u>
UniProt ID:	<u>Q6Y7W8</u>
Synonyms:	2610016F01Rik; A830080H02Rik; A1852361; AW259676; BC006835; mKIAA0642; Tnrc15
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

