

Product datasheet for **MP205799**

Gsx2 Mouse qPCR Primer Pair (NM_133256)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14843
Forward Sequence:	TGTCTGCGCCGCCACTACCTA
Reverse Sequence:	CTGGTAAACGCTGTCCTCATCC
ACCN:	<u>NM_133256</u> , <u>NM_133256.1</u> , <u>NM_133256.2</u> , <u>BC119464</u> , <u>BC119465</u> , <u>BC134359</u>
UniProt ID:	<u>P31316</u>
Synonyms:	Gsh-2; Gsh2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

