

Product datasheet for **MP205805**

Gtpbp3 Mouse qPCR Primer Pair (NM_032544)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	70359
Forward Sequence:	GGTGAGGATTGCGTGGAGTTTC
Reverse Sequence:	GCAAATGCCCTCCTGGTGAAC
ACCN:	<u>NM_032544</u> , <u>NM_032544.1</u> , <u>NM_032544.2</u> , <u>NM_032544.3</u> , <u>BC138560</u> , <u>BC055716</u> , <u>BC075712</u> , <u>BC092042</u> , <u>BC145333</u> , <u>NM_032544.4</u>
UniProt ID:	<u>Q923K4</u>
Synonyms:	2410009F13Rik; A1607903
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

