

Product datasheet for **MP205812**

Gapvd1 Mouse qPCR Primer Pair (NM_025709)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66691
Forward Sequence:	GTCAGAGTCTCTGCTTGCCATG
Reverse Sequence:	CACTGCACCTTCCAGGATTGGA
ACCN:	<u>BC060123</u> , <u>NM_025709</u> , <u>NM_025709.1</u> , <u>NM_025709.2</u> , <u>BC031478</u> , <u>BC032088</u> , <u>BC043715</u> , <u>BC048847</u> , <u>BC057164</u> , <u>NM_025709.3</u>
UniProt ID:	<u>Q6PAR5</u>
Synonyms:	2010005B09Rik; 4432404J10Rik; AW108497; Gapex-5; mKIAA1521; RAP6; RME-6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

