

Product datasheet for **MP205867**

Gnl3l Mouse qPCR Primer Pair (NM_198110)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	237107
Forward Sequence:	CTGCCAACTGTGGCTTTCAAGG
Reverse Sequence:	GAGATTTTCGGCTCCAAAGCAGG
ACCN:	BC057033 , BC079653 , NM_198110 , NM_198110.1 , NM_198110.2 , BB644039 , BC020354 , BY640210
UniProt ID:	Q6PGG6
Synonyms:	BC020354
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)