

## Product datasheet for **MP205868**

### Gdi1 Mouse qPCR Primer Pair (NM\_010273)

#### Product data:

|                   |                                                                                                                                                                                                                                                                                                                                                                   |
|-------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                                                 |
| Gene ID:          | 14567                                                                                                                                                                                                                                                                                                                                                             |
| Forward Sequence: | CCTCAGAACCAAGTCAACAGGA                                                                                                                                                                                                                                                                                                                                            |
| Reverse Sequence: | CTCCAATGCAGGCTCAACTTC                                                                                                                                                                                                                                                                                                                                             |
| ACCN:             | <u><a href="#">BC013758</a></u> , <u><a href="#">BC037598</a></u> , <u><a href="#">NM_010273</a></u> , <u><a href="#">NM_010273.1</a></u> , <u><a href="#">NM_010273.2</a></u> , <u><a href="#">NM_010273.3</a></u> , <u><a href="#">NM_010273.4</a></u> ,<br><u><a href="#">BB613634</a></u> , <u><a href="#">BC010220</a></u> , <u><a href="#">BC025606</a></u> |
| UniProt ID:       | <u><a href="#">P50396</a></u>                                                                                                                                                                                                                                                                                                                                     |
| Synonyms:         | GDIA; GD1alpha                                                                                                                                                                                                                                                                                                                                                    |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.                                                                                                                                                                    |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.                     |
| Storage:          | Store at -20°C.                                                                                                                                                                                                                                                                                                                                                   |
| Stability:        | The primer mix is stable for one year from date of shipping.                                                                                                                                                                                                                                                                                                      |



[View online »](#)