

Product datasheet for **MP205892**

Gulo Mouse qPCR Primer Pair (NM_178747)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	268756
Forward Sequence:	GTTCCACCTTCTGGAGACATCC
Reverse Sequence:	CGTTCTCACTGTGAGGAAACCAG
ACCN:	<u>BC019856</u> , <u>BC028822</u> , <u>BC028828</u> , <u>NM_178747</u> , <u>NM_178747.1</u> , <u>NM_178747.2</u> , <u>NM_178747.3</u> , <u>BC034835</u> , <u>BC035221</u>
UniProt ID:	<u>P58710</u>
Synonyms:	5730581M22; AU018375; BC028822; GLO; LGO; sfx; unh; unhip
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)