

Product datasheet for **MP205894**

Hlf0 Mouse qPCR Primer Pair (NM_008197)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14958
Forward Sequence:	CATGATCGTGGCTGCTATCCAG
Reverse Sequence:	GTTCTCACCCACCTTGTAGTGG
ACCN:	BC003830 , BC011493 , BC110361 , NM_008197 , NM_008197.1 , NM_008197.2 , NM_008197.3 , BB661016
UniProt ID:	P10922
Synonyms:	D130017D06Rik; H1(0); H1-0; H1f; H1fv
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

