

Product datasheet for **MP205904**

H2az1 Mouse qPCR Primer Pair (NM_016750)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51788
Forward Sequence:	GGCCGTATTCATCGACACCTGA
Reverse Sequence:	GACGCATTCCTGCCAACTCAAG
ACCN:	<u>BC079903</u> , <u>NM_016750</u> , <u>NM_016750.1</u> , <u>NM_016750.2</u> , <u>NM_016750.3</u> , <u>BC147478</u> , <u>BC158035</u> , <u>BC171999</u>
UniProt ID:	<u>P0C0S6</u>
Synonyms:	H2A.; H2A.Z; H2a.z-1; H2A.Z1; H2af; H2afz
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)