

Product datasheet for **MP205911**

Hmx3 Mouse qPCR Primer Pair (NM_008257)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	15373
Forward Sequence:	GTCACCTCCCGCGACCAGAAG
Reverse Sequence:	TTGGAGTCCAGCTCCTTGTGGT
ACCN:	<u>NM_008257</u> , <u>NM_008257.1</u> , <u>NM_008257.2</u> , <u>NM_008257.3</u> , <u>BC139482</u> , <u>BC139483</u> , <u>BE649795</u> , <u>BE956569</u>
UniProt ID:	<u>P42581</u>
Synonyms:	Nkx-5.1; Nkx5-1; Nkx5.1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

