

Product datasheet for **MP205950**

Hspa12a Mouse qPCR Primer Pair (NM_175199)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	73442
Forward Sequence:	GCAGAGTCGGACCTTTTGGTG
Reverse Sequence:	CTGGTAACCGTATCTGATGGACG
ACCN:	<u>BC030362</u> , <u>NM_175199</u> , <u>NM_175199.1</u> , <u>NM_175199.2</u> , <u>NM_175199.3</u> , <u>BQ176746</u> , <u>BY240671</u>
UniProt ID:	<u>Q8K0U4</u>
Synonyms:	I700063D12Rik; AI118035; AI840429; AW048913; AW556406; D5Mgi40; Gm19925; mKIAA0417
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

