

Product datasheet for **MP205959**

Hspa4 Mouse qPCR Primer Pair (NM_008300)

Product data:

| | |
|-------------------|---|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 15525 |
| Forward Sequence: | CACGGCATCTTCAGTGTGTCCA |
| Reverse Sequence: | GGTTCTTCTGGTCCACTTGCA |
| ACCN: | <u>BC003770</u> , <u>NM_008300</u> , <u>NM_008300.1</u> , <u>NM_008300.2</u> , <u>NM_008300.3</u> , <u>BC003950</u> |
| Synonyms: | 70kDa; AI317151; APG-2; Hsp70RY; Hsp110; mKIAA4025 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec. |
| Storage: | Store at -20 °C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)